

LifeTec®

Rescue Training Dummies

TRAINING MANIKINS FOR FIRE & RESCUE SERVICES



FIRE & RESCUE SERVICES FOR PROFESSIONALS WHO SAVE LIVES

Ruth Lee manikins are used throughout the world by fire and rescue professionals. They are tough enough for the fire service and military, and are weighted proportionately to be anatomically correct, for more realistic training.

We are confident you won't find a more durable manikin on the market and with a full range of dummies, designed for different training scenarios, we are sure we can provide you with the perfect training partner!

DUTY RANGE

Our best selling manikin and commonly known as the "Drill Dummy". It is ideal for a range of training scenarios including: Confined Space, Urban Search and Rescue (USAR), Road Traffic Collision (RTC), CBRN and general handling exercises.

CODE	DESCRIPTION	WEIGHT	HEIGHT	NATO STOCK NO
RLN5	BABY	5 Kg	0.7m	4240-99-8581987
RLN10	TODDLER	10Kg	0.9m	4240-99-0218800
RLN20	YOUTH	20Kg	1.3m	4240-99-2130513
RLN30	30Kg ADULT	30Kg	1.8m	4240-99-5158249
RLN50	50Kg ADULT	50Kg	1.8m	4240-99-5550594
RLN70	70Kg ADULT	70Kg	1.8m	4240-99-2416796
RLN90	90Kg ADULT	90Kg	1.8m	NA



DUTY PLUS

An enhanced version of our popular Duty manikin. Because of the extra rigidity built into the neck, chest and torso area this manikin acts more like a 'conscious casualty'. It is ideal for vehicle extrication training.

CODE	DESCRIPTION	WEIGHT	HEIGHT	NATO STOCK NO
RLN50+	50 KG ADULT	50Kg	1.8m	NA



CONFINED SPACE

The perfect training partner for confined space training, or any other difficult access location as this manikin splits into several sections making it easier to transport and construct once in position.

CODE	DESCRIPTION	WEIGHT	HEIGHT	NATO STOCK NO
RLNCS	50 KG ADULT	50Kg	1.8m	NA



FIRE HOUSE

This manikin is ideal for snatch rescue training and 'hot cell' training. It will withstand temperatures of around 160C at floor level (400C at ceiling height).

CODE	DESCRIPTION	WEIGHT	HEIGHT	NATO STOCK NO
RLNFH10	TODDLER	10Kg	0.9m	NA
RLNFH20	YOUTH	20Kg	1.3m	NA
RLNFH30	30Kg ADULT	30Kg	1.8m	4240-99-5158249
RLNFH50	50Kg ADULT	50Kg	1.8m	4240-99-5550594
RLNFH70	70Kg ADULT	70Kg	1.8m	NA



FULL BODIED AIRWAY MANAGEMENT

Your team will be expected to carry our airway management skills in some difficult circumstances - not a warm, well-lit classroom! Introducing this manikin into your routine training will allow you to create more realistic scenarios and ensure your team are resilient.


CODE	DESCRIPTION	WEIGHT	HEIGHT
RLNAWM	20KG ADULT WITHOUT TORSO	13Kg	1.8m
RLNAWMTOR	20KG ADULT WITH TORSO	20Kg	1.8m
RLNAWM50	50KG ADULT WITHOUT TORSO	43Kg	1.8m
RLNAWMTOR50	50KG ADULT WITH TORSO	50Kg	1.8m



MULTI TRAUMA

One of our more interesting manikins that can be used to simulate impalements or a partial or total loss of limb. It is an excellent choice for USAR or RTC (Road Traffic Collision) training and helps to promote proper casualty care -after all, if you handle him incorrectly, he could lose an arm!


CODE	DESCRIPTION	WEIGHT	HEIGHT	NATO STOCK NO
RLN30MT	30Kg ADULT	30Kg	1.8m	4210-99-2547551
RLN50MT	50Kg ADULT	50Kg	1.8m	NA
RLN70MT	70Kg ADULT	70Kg	1.8m	NA



WORKING AT HEIGHT

Specially designed for Working at Height training. This manikin is constructed to provide extra rigidity and reduce slumping when using a rescue harness.

CODE	DESCRIPTION	WEIGHT	HEIGHT	NATO STOCK NO
RLWAH30	30Kg ADULT	30Kg	1.8m	NA
RLWAH50	50Kg ADULT	50Kg	1.8m	NA
RLWAH70	70Kg ADULT	70Kg	1.8m	NA



BARIATRIC

Obesity is a growing problem in the Western world. Would you and your team be confident in extricating a heavier casualty in an emergency? With a 'light-weight' option available when bulk and not excessive weight is required.

CODE	DESCRIPTION	WEIGHT	HEIGHT	NATO STOCK NO
RLN90B	BARI ADULT	30Kg	1.8m	NA
RLN180	BADULT	50Kg	1.8m	NA
RLN260	70Kg ADULT	70Kg	1.8m	NA



WATER-FILLABLE BARIATRIC SUIT

A NEW innovation for Bariatric training. Our Water-fillable Bariatric Suit fits over any Duty Adult manikin and allows you to add up to 150kgs of additional weight. When empty it is light enough for one person to carry easily.

CODE	DESCRIPTION
RLBCS	WATER-FILLABLE BARIATRIC SUIT



LifeTec[®]
Rescue Training Dummies

LifeTec Pty Ltd
Unit 44, 6 Jubilee Avenue (J6)
Warriewood
NSW 2102
02 9979 7079
info@lifetec.com.au
lifetec.com.au

